


## Celiac Testing at Prometheus

Order Name: **PRO CELIAC**  
Test Number: 3613275  
Revision Date: 12/08/2016

TEST NAME	METHODOLOGY	LOINC CODE
Celiac Testing at Prometheus		

SPECIMEN REQUIREMENTS				
Specimen	Specimen Volume (min)	Specimen Type	Specimen Container	Transport Environment
Preferred	See Instructions	EDTA Whole Blood and Serum	EDTA (lavender top) and Clot Activator SST	Refrigerated
Instructions	Collect both serum and whole blood. 2mL serum (red/grey or red no-gel) and 5mL whole blood (EDTA). Indicate desired test code: <ul style="list-style-type: none"><li>6360 - Celiac Plus</li><li>6260 - Celiac Genetics HLA DQ2/DQB</li><li>1155 - Celiac Serology</li></ul> Specimen Stability: Room temperature: 7 days, Refrigerated: 30 days, Frozen: Unacceptable			

GENERAL INFORMATION	
Expected TAT	3 days
Notes	Reference Lab: Prometheus
Service Provided By	 <b>labcorp</b> Oklahoma, Inc.