


**Meclizine Screen S/P**

Order Name: **Meclizine Scr**  
Test Number: 6900723  
Revision Date: 03/14/2017

TEST NAME	METHODOLOGY	LOINC CODE
Meclizine Screen S/P	<u>High Performance Liquid Chromatography</u>	

SPECIMEN REQUIREMENTS				
Specimen	Specimen Volume (mL)	Specimen Type	Specimen Container	Transport Environment
Preferred	2.0 mL (0.8 mL)	Plasma	EDTA (Lavender Top)	Refrigerated
Alternate 1	2.0 mL (0.8 mL)	Serum	Clot Activator (Red Top, No-Gel)	Refrigerated
Instructions	Specimen Stability: Room temperature: 14 days, Refrigerated: 1 month, Frozen: 6 months Unacceptable: Collection with gel barriers (SST or PST).			

GENERAL INFORMATION	
Testing Schedule	Tuesday, Friday
Expected TAT	3 days
Clinical Use	Also known as: Antivert, Bonine, RuVert
Notes	Positive screen results in confirmation reflex. Meclizine Confirmation - (Cpt 80375) Order Code: <b>Meclizine Confirmation S/P [6900725]</b> Result Code: <b>Meclizine Con [6900725]</b>  Reference Lab: NMS NMS Test Code: 9196SP <a href="#">Click Here</a> to view information on the NMS website.
CPT Code(s)	80307
Service Provided By	 <b>labcorp</b> Oklahoma, Inc.