


Meclizine Screen S/P

Order Name: **Meclizine Scr**
 Test Number: 6900723
 Revision Date: 03/14/2017

TEST NAME	METHODOLOGY	LOINC CODE
Meclizine Screen S/P	High Performance Liquid Chromatography	

SPECIMEN REQUIREMENTS				
Specimen	Specimen Volume (min)	Specimen Type	Specimen Container	Transport Environment
Preferred	2.0 mL (0.8 mL)	Plasma	EDTA (Lavender Top)	Refrigerated
Alternate 1	2.0 mL (0.8 mL)	Serum	Clot Activator (Red Top, No-Gel)	Refrigerated
Instructions	Specimen Stability: Room temperature: 14 days, Refrigerated: 1 month, Frozen: 6 months Unacceptable: Collection with gel barriers (SST or PST).			

GENERAL INFORMATION	
Testing Schedule	Tuesday, Friday
Expected TAT	3 days
Clinical Use	Also known as: Antivert, Bonine, RuVert
Notes	Positive screen results in confirmation reflex. Meclizine Confirmation - (Cpt 80375) Order Code: Meclizine Confirmation S/P [6900725] Result Code: Meclizine Con [6900725] Reference Lab: NMS NMS Test Code: 9196SP Click Here to view information on the NMS website.
CPT Code(s)	80307
Service Provided By	 Oklahoma, Inc.