

**EBV (Epstein Barr Virus) Viral Capsid Antibodies IgG and IgM**

Order Name: **EBV VCA Ab G/M**  
 Test Number: 5580926  
 Revision Date: 04/11/2023

TEST NAME	METHODOLOGY	LOINC CODE
<b>EBV VCA IgG</b>	<u>Chemiluminescence Assay</u>	5157-3
<b>EBV VCA IgG Interpretation</b>		24114-1
<b>EBV VCA IgM</b>	<u>Chemiluminescence Assay</u>	5159-9
<b>EBV VCA IgM Interpretation</b>		24115-8

**SPECIMEN REQUIREMENTS**

Specimen	Specimen Volume (min)	Specimen Type	Specimen Container	Transport Environment
Preferred	<b>2 mL (0.25mL)</b>	<b>Serum</b>	<b>Clot Activator SST</b>	<b>Refrigerated or Frozen</b>
Instructions	Specimen Stability: Room Temperature: Unacceptable, Refrigerated: 7 days, Frozen 30 Days			

**GENERAL INFORMATION**

Testing Schedule	Mon-Sat
Expected TAT	1-2 Days
Notes	EB Virus IgM IFA confirmation testing will no longer be performed due to supply chain availability. The results for this test are not, by themselves, diagnostic and should be considered in association with clinical presentation and clinical history.
CPT Code(s)	86665x2
Service Provided By	 <b>labcorp</b> Oklahoma, Inc.