


Herpes Select Panel (HSV 1/2 IgG and IgM)

Order Name: **Herpes Panel**
Test Number: **5563986**
Revision Date: **04/20/2020**

TEST NAME	METHODOLOGY	LOINC CODE
HSV Type 1 IgG	<u>Chemiluminescence Assay</u>	5206-8
HSV Type 1 IgG Interpretation		51916-5
HSV Type 2 IgG	<u>Chemiluminescence Assay</u>	5209-2
HSV Type 2 IgG Interpretation		43180-9
<u>Herpes Simplex Virus Type 1/2 IgM IFA</u>	<u>Indirect Fluorescent Antibody</u>	

SPECIMEN REQUIREMENTS				
Specimen	Specimen Volume (min)	Specimen Type	Specimen Container	Transport Environment
Preferred	3 mL	Serum	Clot Activator SST	Refrigerated or Frozen
Instructions	Specimen Stability: Room Temperature: Unacceptable, Refrigerated: 7 days, Frozen 30 Days			

GENERAL INFORMATION	
Testing Schedule	Varies
Expected TAT	3 - 5 days
Notes	HSV 1/2 IgG are performed Mon - Sat and HSV 1/2 IgM are performed Tues/Fri
CPT Code(s)	86695x2, 86696x2
Service Provided By	 labcorp Oklahoma, Inc.