


Deoxypyridinoline Crosslinks (Pyrilinks-D)

Order Name: **Pyrilinks-D**
Test Number: 6904991
Revision Date: 01/28/2025

TEST NAME	METHODOLOGY	LOINC CODE
Deoxypyridinoline Crosslinks (Pyrilinks-D)	<u>Immunochemiluminometric (ICMA)</u>	

SPECIMEN REQUIREMENTS				
Specimen	Specimen Volume (mL)	Specimen Type	Specimen Container	Transport Environment
Preferred	10 mL (1 mL)	Urine, First Void	Sterile Urine container	Refrigerated
Instructions	<p>Specimen Type: Plastic urine container Specimen Storage: Refrigerated Specimen Collection: Collection before 10 AM is recommended (eg, first morning void). Do NOT add preservative. Special Instructions: The account must submit the patient's Social Security number to monitor. Specimen Stability: Ambient: 14 days, Refrigerated : 14 days, Frozen: 14 days</p>			

GENERAL INFORMATION	
Expected TAT	4-5 days
CPT Code(s)	82570, 82523
Service Provided By	 labcorp Oklahoma, Inc.