


**HLA-B\*5701 Typing (Abacavir hypersensitivity)**

Order Name: **HLA-B5701**  
Test Number: **6905039**  
Revision Date: **12/09/2022**

| TEST NAME                                     | METHODOLOGY                      | LOINC CODE |
|---|----------------------------------|------------|
| HLA-B*5701 Typing (Abacavir hypersensitivity) | <u>Polymerase Chain Reaction</u> |            |

| SPECIMEN REQUIREMENTS |  |               |                     |                       |
|-----------------------|--|---------------|---------------------|-----------------------|
| Specimen              | Specimen Volume (mL)   | Specimen Type | Specimen Container  | Transport Environment |
| Preferred             | 7 mL (3 mL)  | Whole Blood   | EDTA (Lavender Top) | Room Temperature      |
| Instructions          | <p><b>Notes:</b> 3 mL whole blood or four (4) buccal swabs (Note: This volume Does NOT allow for repeat testing.)</p> <p><b>Specimen Type:</b> Lavender-top (EDTA) tube or four (4) buccal swabs in a sealed envelope (buccal swab kit). When submitting buccal swabs please use a buccal swab kit provided by LabCorp. To obtain the buccal kit, please telephone 800-533-1037.</p> <p><b>Container Detail:</b> WhiBld</p> <p><b>Specimen Storage:</b> Maintain blood specimen at room temperature. Keep buccal swabs dry and at room temperature.</p> <p><b>Specimen Collection:</b> Not Available</p> <p><b>Special Instructions:</b> If you have questions, please telephone 800-533-1037 (HLA customer service) for assistance selecting the proper HLA test for the patient.</p> <p><b>Specimen Stability:</b> Ambient: Not Available, Refrigerated : Not Available, Frozen: Not Available</p> |               |                     |                       |

| GENERAL INFORMATION |  |
|---------------------|--|
| Expected TAT        | 3 - 7 days   |
| Notes               | Labcorp Test Code: 006926  |
| CPT Code(s)         | 81381  |
| Service Provided By |  <b>labcorp</b><br>Oklahoma, Inc. |