

Infliximab Activity and Neutralizing Antibody

Order Name: **Infliximab NAB**
Test Number: 6905051
Revision Date: 01/02/2020

| TEST NAME | METHODOLOGY | LOINC CODE |
|---|-------------------------------------|------------|
| Infliximab Activity and Neutralizing Antibody | Electrochemiluminescent Immunoassay | |

| SPECIMEN REQUIREMENTS | | | | |
|-----------------------|--|---------------|--------------------|-----------------------|
| Specimen | Specimen Volume (min) | Specimen Type | Specimen Container | Transport Environment |
| Preferred | 2 mL (1 mL) | Serum | Clot Activator SST | Frozen |
| Instructions | <p>Notes: 1 mL (Note: This volume Does NOT allow for repeat testing).</p> <p>Specimen Type: Red-top tube or gel-barrier tube</p> <p>Specimen Storage: Frozen (preferred) or refrigerated Stability: Room temperature: 7 days Refrigerated: 7 days Frozen: 359 days Freeze/thaw cycle: stable x6</p> <p>Specimen Collection: Serum must be separated from cells within 45 minutes of venipuncture. Send serum in a plastic transport tube. To avoid delays in turnaround time when requesting multiple test on frozen samples, PLEASE SUBMIT SEPARATE FROZEN SPECIMENS FOR EACH TEST REQUESTED.</p> <p>Special Instructions: Testing is performed at Esoterix Endocrinology Laboratory: ES#803870 Test includes Serial Monitoring.</p> | | | |

| GENERAL INFORMATION | |
|---------------------|--|
| Expected TAT | 5-10 days |
| Notes | Labcorp Test Code: 503870 |
| CPT Code(s) | 80230, 82397 |
| Service Provided By |  labcorp Oklahoma, Inc. |