

**Trichomonas PAP - TMA**

Order Name: **Trich PAP**  
Test Number: **6910119**  
Revision Date: **10/10/2018**

TEST NAME	METHODOLOGY	LOINC CODE
Trichomonas PAP	Transcription-Mediated Amplification	70166-4

SPECIMEN REQUIREMENTS				
Specimen	Specimen Volume (min)	Specimen Type	Specimen Container	Transport Environment
Preferred	See Instructions	PAP specimen	ThinPrep PreservCyt solution	Room Temperature
Instructions	Cervical brush/spatula in ThinPrep Pap test collection kit. ThinPrep PreservCyt solution. <b>Stability:</b> Specimen must be transferred from ThinPrep PAP vial into an Aptima Specimen Transfer tube within 30 days. Once transferred, the specimen is stable in the Aptima vial for 30 days refrigerated.  <b>Please contact your account rep if training and ThinPrep supplies have not been supplied to your practice.</b>			

GENERAL INFORMATION	
Testing Schedule	Mon-Fri
Expected TAT	1-2 days after set-up
CPT Code(s)	87661
Service Provided By	 <b>labcorp</b> Oklahoma, Inc.