


Donath-Landsteiner Test

Order Name: **Donath-Land**
Test Number: **7306200**
Revision Date: **09/24/2019**

TEST NAME	METHODOLOGY	LOINC CODE
Donath-Landsteiner Test	Hemolysis	

SPECIMEN REQUIREMENTS				
Specimen	Specimen Volume (min)	Specimen Type	Specimen Container	Transport Environment
Preferred	10 mL (7 mL)	Serum	Clot Activator (Red Top, No-Gel)	Refrigerated
Instructions	Blood should be allowed to clot at 37°C, and the serum separated from the clot at this temperature. Once the serum is separated from the cells, the serum can be stored and transported refrigerated.			

GENERAL INFORMATION	
Testing Schedule	Mon-Fri (7:00am - 3:30pm)
Expected TAT	1 Day
Clinical Use	Used in the diagnosis of Paroxysmal Cold Hemoglobinuria (PCH).
CPT Code(s)	86940 86941
Service Provided By	 labcorp Oklahoma, Inc.