


Pepper C. annum (F218-IgE Paprika/Sweet Pepper)

Order Name: **Pepper C Annum**
Test Number: **5613535**
Revision Date: **03/24/2025**

TEST NAME	METHODOLOGY	LOINC CODE
Pepper C. annum (F218-IgE Paprika/Sweet Pepper)	ImmunoCAP	

SPECIMEN REQUIREMENTS				
Specimen	Specimen Volume (min)	Specimen Type	Specimen Container	Transport Environment
Preferred	0.25 mL (0.1 mL)	Serum	Clot Activator SST	Room Temperature
Instructions	Specimen Type: Red-top tube or gel-barrier tube Specimen Storage: Room temperature Specimen Collection: Not Available Specimen Stability: Ambient: 14 days, Refrigerated : 14 days, Frozen: 3 months			

GENERAL INFORMATION	
Expected TAT	2 - 4 days
Clinical Use	The term sweet pepper covers a wide variety of mild peppers that, like the chile, belong to the Capsicum family. The best known sweet peppers are bell peppers. This family also includes paprika, green pepper, Hungarian pepper, red pepper and pimento. Please note that this allergen is not related to the peppercorn spices (green and black pepper).
Notes	Labcorp Test Code: 602779
CPT Code(s)	86003
Service Provided By	 labcorp Oklahoma, Inc.