

Yellow Dock IgE

Order Name: **Yellow Dock**
Test Number: **5613525**
Revision Date: **12/14/2022**

TEST NAME	METHODOLOGY	LOINC CODE
Yellow Dock IgE	ImmunoCAP	

SPECIMEN REQUIREMENTS				
Specimen	Specimen Volume (min)	Specimen Type	Specimen Container	Transport Environment
Preferred	0.25 mL (0.1 mL)	Serum	Clot Activator SST	Room Temperature
Instructions	Specimen Type: Red-top tube or gel-barrier tube Specimen Storage: Room temperature Specimen Collection: Not Available Specimen Stability: Ambient: 14 days, Refrigerated : 14 days, Frozen: 3 months			

GENERAL INFORMATION	
Expected TAT	3 - 4 days
Notes	Labcorp Test Code: 602974
CPT Code(s)	86003
Service Provided By	 labcorp Oklahoma, Inc.