

Benzodiazepines, Serum or Plasma, Quantitative

Order Name: **Benzodiaz S/P**
 Test Number: 6980151
 Revision Date: 01/04/2023

TEST NAME	METHODOLOGY	LOINC CODE
Benzodiazepines, Serum or Plasma, Quantitative	Quantitative Liquid Chromatography/Tandem Mass Spectrometry	

SPECIMEN REQUIREMENTS

Specimen	Specimen Volume (min)	Specimen Type	Specimen Container	Transport Environment
Preferred	3 mL (1 mL)	Serum	Clot Activator (Red Top, No-Gel)	Refrigerated
Instructions	Specimen Type: Red-top tube OR green-top (heparin) tube Specimen Storage: SUBMISSION/TRANSPORT (<3 days): Room Temperature. For storage beyond three days, specimen should be refrigerated or frozen. Specimen Collection: Serum or plasma should be separated from cells within two hours of venipuncture. Send serum or plasma in a plastic transport tube. Special Instructions: Testing referred to MEDTOX Laboratories Inc TC19976			

GENERAL INFORMATION

Expected TAT	5-10 days
Notes	Labcorp Test Code: 791886
CPT Code(s)	80347 (Alt code: G0480)
Service Provided By	 labcorp Oklahoma, Inc.