

## MySwab Vaginosis Profile + Trichomonas - Swab - TMA

Order Name: **MySwab Vag TV**  
 Test Number: 6987004  
 Revision Date: 04/06/2020

TEST NAME	METHODOLOGY	LOINC CODE
Bacterial Vaginosis - Swab	<a href="#">Transcription-Mediated Amplification</a>	46999-9
Candida Species - Swab	<a href="#">Transcription-Mediated Amplification</a>	47000-5
Candida Glabrata - Swab	<a href="#">Transcription-Mediated Amplification</a>	n/a
Trichomonas Swab (Vaginal)	<a href="#">Transcription-Mediated Amplification</a>	70165-6

SPECIMEN REQUIREMENTS				
Specimen	Specimen Volume (min)	Specimen Type	Specimen Container	Transport Environment
Preferred	<b>Swab</b>	<b>Swab</b>	<b>APTIMA Multitest Swab kit</b>	<b>Refrigerated</b>
Instructions	Collect Swab in APTIMA Multitest Swab collection kit. Place swab in Swab Specimen Transport Tube, break shaft off at scoreline then recap tube. Stability: Room Temperature 30 Days, Refrigerated 30 Days, Frozen: N/A			

GENERAL INFORMATION	
Testing Schedule	Mon-Fri
Expected TAT	1-2 days after set-up
CPT Code(s)	87801, 87481, 87661
Service Provided By	 Oklahoma, Inc.