

**Bacterial Vaginosis - Swab - TMA**

Order Name: **Bact Vag Swab**  
Test Number: 6987001  
Revision Date: 01/01/2021

TEST NAME	METHODOLOGY	LOINC CODE
Bacterial Vaginosis - Swab	Transcription-Mediated Amplification	46999-9

SPECIMEN REQUIREMENTS				
Specimen	Specimen Volume (min)	Specimen Type	Specimen Container	Transport Environment
Preferred	Swab	Swab	APTIMA Multitest Swab kit	Ambient / Refrigerated
Instructions	Collect Swab in APTIMA Multitest Swab collection kit. Place swab in Swab Specimen Transport Tube, break shaft off at scoreline then recap tube. Stability: Room Temperature 30 Days, Refrigerated 30 Days, Frozen: N/A			

GENERAL INFORMATION	
Testing Schedule	Mon-Fri
Expected TAT	1-2 days after set-up
CPT Code(s)	81513
Service Provided By	 <b>labcorp</b> Oklahoma, Inc.