


Peripheral Blood for Mycosis Fungoides by Flow

Order Name: **Mycosis Flow**
Test Number: 5582601
Revision Date: 07/06/2022

TEST NAME	METHODOLOGY	LOINC CODE
Peripheral Blood for Mycosis Fungoides by Flow	Flow cytometry	

SPECIMEN REQUIREMENTS				
Specimen	Specimen Volume (mL)	Specimen Type	Specimen Container	Transport Environment
Preferred	4.5 mL (1 mL)	Whole Blood	EDTA (Lavender Top)	Room Temperature
Instructions	Deliver to laboratory (flow cytometry) ASAP. DO NOT Centrifuge or Refrigerate. Collect Monday through Friday only! Specimen must be received in the lab section by Friday afternoon or collect sample on Monday. Specimen stability: 48hrs Room Temperature.			

GENERAL INFORMATION	
Testing Schedule	Tue-Sat
Expected TAT	2 days
Notes	Results from this test will be reported in a pathology report: Order Code: [7120723] Flow Cytometry Final Report Result Code: [7120723] Flow Cytometry Final Report
CPT Code(s)	88184, 88185X18, 88187
Service Provided By	 labcorp Oklahoma, Inc.