


Anti-Hu, Anti-Ri, Anti-Yo IFA

Order Name: **Hu Ri Yo IFA**
Test Number: **5194477**
Revision Date: **09/27/2023**

TEST NAME	METHODOLOGY	LOINC CODE
Anti-Hu, Anti-Ri, Anti-Yo IFA	Immunofluorescence Assay (IFA)	

SPECIMEN REQUIREMENTS				
Specimen	Specimen Volume (min)	Specimen Type	Specimen Container	Transport Environment
Preferred	2 mL (1 mL)	Serum	Clot Activator (SST or Red No-Gel)	Room Temperature
Instructions	<p>Notes: 1 mL (Note: This volume Does not allow for repeat testing.)</p> <p>Specimen Type: Serum from a Red Top Tube, Gel-barrier Tube, Serum Transfer tube</p> <p>Specimen Storage: Room Temperature</p> <p>Specimen Collection: Not Available</p> <p>Special Instructions: If reflex test is performed, additional charges/CPT code(s) may apply.</p> <p>Specimen Stability: Ambient: 14 Days, Refrigerated : 14 Days, Frozen: 14 Days</p>			

GENERAL INFORMATION	
Expected TAT	3 - 5 days
CPT Code(s)	86255x3 If positive a reflex test confirmatormation is performed, additional charges/CPT code(s) may apply. .Rflx Hu Blot [5194481] (.Reflex Anti-Hu Ab Line Blot Rflx Result code [5194481]) cpt: 84182 .Rflx Ri Blot [5194483] (.Reflex Anti-Ri Ab Line Blot Rflx Result code [5194483]) cpt: 84182 .Rflx Yo Blot [5194485] (.Reflex Anti-Yo Ab Line Blot Rflx Result code [5194485]) cpt: 84182
Service Provided By	 labcorp Oklahoma, Inc.