



Parainfluenza 1,2,3 by Complement Fixation

Order Name: **PARAINFLU**
Test Number: 5504941
Revision Date: 12/15/2022

TEST NAME	METHODOLOGY	LOINC CODE
Parainfluenza 1,2,3 by Complement Fixation	Complement fixation (CF)	1668-3

SPECIMEN REQUIREMENTS

Specimen	Specimen Volume (min)	Specimen Type	Specimen Container	Transport Environment
Preferred	2 mL	Serum	Clot Activator SST	Refrigerated
Instructions	Specimen Type: Red-top tube or gel-barrier tube Specimen Storage: Refrigerate Specimen Collection: Not Available Specimen Stability: Ambient: Not Available, Refrigerated : Not Available, Frozen: Not Available			

GENERAL INFORMATION

Expected TAT	1 - 5 days 2 days.
CPT Code(s)	86710x3
Service Provided By	 labcorp Oklahoma, Inc.