


**Loxapine Level**

Order Name: **LOXAPINE**  
Test Number: 4001200  
Revision Date: 11/06/2025

TEST NAME	METHODOLOGY	LOINC CODE
Loxapine Level	<u>Liquid Chromatography/Tandem Mass Spectrometry</u>	

SPECIMEN REQUIREMENTS				
Specimen	Specimen Volume (min)	Specimen Type	Specimen Container	Transport Environment
Preferred	1 mL (0.3 mL)	Serum	Clot Activator (Red Top, No-Gel)	Refrigerated
Alternate 1	1 mL (0.3 mL)	Plasma	Lithium Heparin (Dark Green Top / No-Gel)	Refrigerated
Instructions	<p><b>Specimen:</b> Gel-barrier tubes are NOT acceptable.</p> <p><b>Volume:</b> 1 mL (0.3 mL) Serum Red-Top Tube or Plasma Lavender-Top (EDTA) Tube</p> <p><b>Collection:</b> Promptly centrifuge and separate Serum or Plasma into a plastic screw capped vial using approved guidelines. Testing referred to National Med Services NMSPA#2538SP</p> <p><b>Stability Requirements:</b> Room temperature 30 days, Refrigerated 30 days, Frozen 15 months</p> <p>Cause for Rejection: Gel Barrier Tube</p>			

GENERAL INFORMATION	
Expected TAT	11-13 days
CPT Code(s)	80342
Service Provided By	 <b>labcorp</b> Oklahoma, Inc.