

Ovarian Malignancy Risk-ROMA

Order Name: **Ovarian Risk**
Test Number: **5195163**
Revision Date: **10/30/2024**

TEST NAME	METHODOLOGY	LOINC CODE
CA 125 Assay	Chemiluminescence Assay	10334-1
HE4, Ovarian Cancer Monitoring	Quantitative Electrochemiluminescent Immunoassay	55180-4
Roma Premenopausal ROMA	Interpretive information.	69569-2
Roma Postmenopausal ROMA	Interpretive information.	69570-0
Roma Risk Comment	Interpretive information.	77202-0

SPECIMEN REQUIREMENTS				
Specimen	Specimen Volume (min)	Specimen Type	Specimen Container	Transport Environment
Preferred	1.75 mL (1 mL)	Serum	Clot Activator (SST or Red No-Gel)	Refrigerated
Instructions	Specimen: Serum Storage requirements: Refrigerate; stable for 5 days at room temperature or 14 days refrigerated. ?Causes for Rejection: Gross hemolysis			