



First Trimester Screen With NT

Order Name: **1st Tri Scrn NT**
Test Number: 5195164
Revision Date: 06/02/2023

| TEST NAME | METHODOLOGY | LOINC CODE |
|--------------------------------|--------------------------------|------------|
| First Trimester Screen With NT | Chemiluminescence, Immunoassay | N/A |

SPECIMEN REQUIREMENTS

| | | | | |
|--------------|---|---------------|--------------------|-----------------------|
| Specimen | Specimen Volume (min) | Specimen Type | Specimen Container | Transport Environment |
| Preferred | 3mL | Serum | Clot Activator SST | Room Temperature |
| Instructions | Specimen: Serum Container: Gel-barrier tube Storage requirements: Room Temperature | | | |

GENERAL INFORMATION

| | |
|---------------------|--|
| Expected TAT | 2-5 days |
| CPT Code(s) | 84163; 84702; 86336 |
| Service Provided By |  labcorp Oklahoma, Inc. |