

## ToxAssure Flex 22, Urine

Order Name: **TA Flex 22 UR**  
 Test Number: 5195611  
 Revision Date: 01/09/2024

TEST NAME	METHODOLOGY	LOINC CODE
ToxAssure Flex 22, Urine	<u>Immunoassay (IA)</u>	

### SPECIMEN REQUIREMENTS

Specimen	Specimen Volume (min)	Specimen Type	Specimen Container	Transport Environment
Preferred	<b>30 mL (15 mL)</b>	<b>Urine, Random</b>	<b>Sterile Urine Container Preservative Free</b>	<b>Room Temperature</b>
<b>Instructions</b>	Specimen: 30mL(15mL) Random Urine Preservative free Urine Transport tube Room Temperature or Refrigerated. Cause for Rejection: Insufficient volume, no ID on container, urine from preservative tube. SUBMISSION/TRANSPORT (<3 days): Room Temperature. For storage beyond three days, specimen should be refrigerated or frozen. Stability: Room Temperature 14 days, Refrigerated 14 days, Frozen 6 weeks.			

### GENERAL INFORMATION

Expected TAT	5-8 days
CPT Code(s)	80347, 80357, 80354, 80307, 80362, 80348, 80368
Service Provided By	