

Encephalitis Panel, Serum with Reflex

Order Name: **Enceph Pnl S**
Test Number: 5197022
Revision Date: 11/20/2024

TEST NAME	METHODOLOGY	LOINC CODE
West Nile Virus Ab,IgG	Semi-Quantitative ELISA/Semi-Quantitative CLIA	38997-3
West Nile Virus Ab,IgM	Semi-Quantitative ELISA/Semi-Quantitative CLIA	38166-5
Mumps Virus Ab IgG	Semi-Quantitative ELISA/Semi-Quantitative CLIA	25418-5
Mumps Virus Ab IgM	Semi-Quantitative ELISA/Semi-Quantitative CLIA	25419-3
Varicella-Zoster Virus Ab,IgG	Semi-Quantitative ELISA/Semi-Quantitative CLIA	5403-1
Varicella-Zoster Virus Ab,IgM	Semi-Quantitative ELISA/Semi-Quantitative CLIA	5404-9
Measles,Rubeola,Antibody IgG	Semi-Quantitative ELISA/Semi-Quantitative CLIA	7962-4
Measles,Rubeola,Antibody IgM	Semi-Quantitative ELISA/Semi-Quantitative CLIA	25421-9
HSV Type 1/2 Combined Ab,IgG	Semi-Quantitative ELISA/Semi-Quantitative CLIA	27948-9

SPECIMEN REQUIREMENTS

Specimen	Specimen Volume (min)	Specimen Type	Specimen Container	Transport Environment
Preferred	4 mL (2 mL)	Serum	Clot Activator SST	Refrigerated
Instructions	Specimen: 4mL(2mL) Serum from Gel Barrier Clot Tube - Refrigerated Collection: Transfer serum to standard transport tube Testing referred to ARUP Y8#3017751 Stability Requirements: Room temperature 48 hours, Refrigerated 2 weeks, Frozen 1 month. Cause for Rejection: CSF. Contaminated, heat-inactivated, or hemolyzed specimens			

GENERAL INFORMATION

Expected TAT	3-7 days
Notes	If HSV 1 and/or 2 IgG is 1.10 IV or greater, then HSV 1 G-specific IgG and HSV 2 G-specific IgG will be added. Additional charges apply. Order Code: .Rflx HSV 1 Glycoprotein G Ab,IgG [5197032] (Result code HSV 1 Gly G IgG [5197032]) cpt: 86695 Order Code: .Rflx HSV 2 Glycoprotein G Ab,IgG [5197033] (Result code HSV 2 Gly G IgG [5197033]) cpt: 86696
CPT Code(s)	86765x2; 86735x2; 86787x2; 86789; 86788; 86694; if reflexed, add 86695; 86696
Service Provided By	 labcorp Oklahoma, Inc.