

Eosinophil Count, Absolute

Order Name: **EOS CT ABS**Test Number: 0100050
Revision Date: 10/31/2019

TEST NAME			METHODOLOGY	LOINC CODE	
Eosinophil Count, Absolute			Flow cytometry	711-2	
SPECIMEN REQUIREMENTS					
Specimen	Specimen Volume (min)	Specimen Type	Specimen Container	Transport Environment	
Preferred	5 mL (1 mL)	Whole Blood	EDTA (Lavender Top)	Refrigerated	
Alternate 1	1 mL (0.5 mL)	Whole Blood	EDTA (Lavender) Microtainer/Bullet	Refrigerated	
Instructions	For best results Room temperature specimens Specimens received greater th		otherwise send Refrigerated. Refrigerated specime	ens can be tested up to 48hrs.	

GENERAL INFORMATION	
Testing Schedule	Daily
Expected TAT	1 Day
Clinical Use	Eosinophilia is found in parasitic diseases, allergic diseases, skin disorders, and certain malignancies.
CPT Code(s)	85004
Service Provided By	labcorp Oklahoma, Inc.