

Phosphorus Urine Random

Order Name: PHOS R U
Test Number: 3002300
Revision Date: 05/28/2021

| TEST NAME | | | METHODOLOGY | LOINC CODE | |
|-------------------------|-----------------------|---|-------------------------|-----------------------|--|
| Phosphorus Urine Random | | Phosphomolybdate Complex | 2778-9 | | |
| SPECIMEN REQUIREMENTS | | | | | |
| Specimen | Specimen Volume (min) | Specimen Type | Specimen Container | Transport Environment | |
| Preferred | 1 mL (0.5) | Urine, Random | Sterile Urine container | Refrigerated | |
| Instructions | | Random urine collection. Diurnal variation exists. Specimen stability: Ambient 8 hours without acid (2 days if acidified <5 pH). | | | |

| GENERAL INFORMATION | |
|---------------------|---|
| Testing Schedule | Daily |
| Expected TAT | 1-2 days |
| Clinical Use | Useful in the diagnosis and treatment of various disorders including parathyroid gland and kidney diseases and vitamin D imbalance. |
| CPT Code(s) | 84105 |
| Service Provided By | labcorp Oklahoma, Inc. |