

Serous Fluid

Order Name: SRS AMYLSE Test Number: 3500050 Revision Date: 04/26/2021

| TEST NAME | | | METHODOLOGY | LOINC CODE | |
|-----------------------|-----------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------|----------------------------|-----------------------|--|
| Amylase Serous Fluid | | | CNPG3 | 1795-4 | |
| SPECIMEN REQUIREMENTS | | | | | |
| Specimen | Specimen Volume (min) | Specimen Type | Specimen Container | Transport Environment | |
| Preferred | 1 mL (0.5) | Serous Fluid | Sterile Screwtop Container | Refrigerated | |
| Instructions | | Venous blood is often drawn simultaneously. Note fluid type on requisition and container. Specimen stability: Ambient=7 days; Refrigerated= 7 days | | | |

| GENERAL INFORMATION | | |
|---------------------|---------------------------------------------------------|--|
| Testing Schedule | Daily | |
| Expected TAT | 1-2 days | |
| Clinical Use | Diagnosis and useful for monitoring acute pancreatitis. | |
| CPT Code(s) | 82150 | |
| Service Provided By | Oklahoma, Inc. | |