

Theophylline

Order Name: **THEOPHYLLI**Test Number: 4005050
Revision Date: 06/04/2021

TEST NAME		MET	THODOLOGY	LOINC CODE
Theophylline		Enz	zyme Immunoassay	4049-3
SPECIMEN REQUIREM	MENTS			
Specimen	Specimen Volume (min)	Specimen Type	Specimen Container	Transport Environment
Preferred	1 mL (0.5)	Plasma	Lithium Heparin PST (Light Green Top)	Refrigerated
Alternate 1	1 mL (0.5)	Serum	Clot Activator SST	Refrigerated
Instructions	Peak: Draw 2-4 hours after oral dose. Stability: Ambient 8 hours, Refrigerated 3 months, Frozen 3 months.			

GENERAL INFORMATION		
Testing Schedule	Daily	
Expected TAT	1-2 days	
Clinical Use	Useful for optimizing drug dosage and assessing toxicity.	
CPT Code(s)	80198	
Service Provided By	labcorp Oklahoma, Inc.	