


Prograf (FK506) Tacrolimus

Order Name: **PROGRAF**
Test Number: 4503275
Revision Date: 05/28/2021

TEST NAME	METHODOLOGY	LOINC CODE
Prograf (FK506) Tacrolimus	Chemiluminescent Micropartical Immunoassay	11253-2

SPECIMEN REQUIREMENTS				
Specimen	Specimen Volume (min)	Specimen Type	Specimen Container	Transport Environment
Preferred	1 mL (0.5)	Whole Blood	EDTA (Lavender Top)	Refrigerated
Instructions	To be drawn 12 hours after dose for 12 hour trough or 24 hour after dose for 24 hour trough. Whole Blood Stability: Stability: Ambient 1 day. Refrigerated 7 days, Frozen 6 months			

GENERAL INFORMATION	
Testing Schedule	Daily
Expected TAT	24hrs (same day if specimen is in lab by 11am)
Clinical Use	Useful for assessing the adequacy of systemic drug delivery since metabolism can exhibit significant variability.
Notes	Also known as Tacrolimus
CPT Code(s)	80197
Service Provided By	 labcorp Oklahoma, Inc.