

**Sodium, Feces**

Order Name: **SODIUM FEC**  
Test Number: 3503100  
Revision Date: 12/12/2022

TEST NAME	METHODOLOGY	LOINC CODE
Sodium, Feces	<u>Ion-Selective Electrode</u>	

SPECIMEN REQUIREMENTS				
Specimen	Specimen Volume (min)	Specimen Type	Specimen Container	Transport Environment
Preferred	5g (1g)	Fecal/Stool	Sterile Orange Screwtop Container	Refrigerated
Instructions	<b>Specimen Type:</b> Stool container (Para Pak-White) <b>Specimen Storage:</b> Refrigerated <b>Specimen Collection:</b> Do not add saline or water to liquefy specimen. <b>Special Instructions:</b> Testing referred to ARUP TC 20379 <b>Specimen Stability:</b> Ambient: 1 hour, Refrigerated : 2 weeks, Frozen: 6 months			

GENERAL INFORMATION	
Expected TAT	5 - 9 days
Notes	Labcorp Test Code: 820607
CPT Code(s)	84302
Service Provided By	 <b>labcorp</b> Oklahoma, Inc.