

Kidney Stone Analysis

Order Name: **STONE ANY**
Test Number: 9101850
Revision Date: 11/15/2021

TEST NAME	METHODOLOGY	LOINC CODE
Kidney Stone Analysis	<u>Infrared Spectroscopy</u>	

SPECIMEN REQUIREMENTS				
Specimen	Specimen Volume (min)	Specimen Type	Specimen Container	Transport Environment
Preferred	n/a	Stone	Sterile Screwtop Container	Room Temperature
Instructions	<p>Specimen Type: 30/50mL clean specimen container or capsule</p> <p>Specimen Storage: Maintain specimen at room temperature</p> <p>Specimen Collection: Provide source (Urinary, gall, salivary) of calculi. Avoid rejection and delayed processing by sending the stone completely dried in a clean container. Do not send gauze, filters, tape or Q-tips. Stones that are not received completely dried will result in delayed processing.</p> <p>Special Instructions: Testing performed at Litholink Stone Analysis Laboratory</p>			

GENERAL INFORMATION	
Expected TAT	4-6 days
Notes	Labcorp Test Code: 910180
CPT Code(s)	82365
Service Provided By	 labcorp Oklahoma, Inc.