

**Estrogens Fractionated** 

Order Name: ESTROG FRA

Test Number: 3618875 Revision Date: 02/22/2023

TEST NAME			METHODOLOGY	LOINC CODE	
Estrogens Fractionated			Liquid Chromatography/Mass Spectroscopy		
SPECIMEN REQU	JIREMENTS				
Specimen	Specimen Volume (min)	Specimen Type	Specimen Container	Transport Environment	
Preferred	2.6 mL (1.3 mL)	Serum	Clot Activator SST	Refrigerated	

Instructions	Gel-barrier serum separator tube. Red top tubes (without gel barrier) should not be used.  1.3 mL (NOTE:This volume does NOT allow for repeat testing.)
<b>GENERAL INFORMAT</b>	ION

GENERAL INFORMATION		
Expected TAT	Report available: 1-4 days from set up at performing lab	
Clinical Use	Estriol Level Estradiol Ultra Sensitive Estrone Level	
Notes	Labcorp Test Code: 146068	
CPT Code(s)	82677, 82679, 82670	
Service Provided By	Oklahoma, Inc.	