

Pediatric Environmental Allergen Panel

Order Name: **A PED EN P**
Test Number: 5622600
Revision Date: 09/11/2023

TEST NAME	METHODOLOGY	LOINC CODE
Aspergillus fumigatus Allergen	ImmunoCAP	15549-9
Cat Dander Allergen	ImmunoCAP	15609-1
Dog Dander Allergen	ImmunoCAP	15685-1
Alternaria alternata Allergen	ImmunoCAP	15530-9
Cladosporium herbarum Allergen	ImmunoCAP	6075-6
House Dust Mite Allergen	ImmunoCAP	15682-8
Dust Mite (D. farinae) Allergen	ImmunoCAP	15680-2

SPECIMEN REQUIREMENTS				
Specimen	Specimen Volume (min)	Specimen Type	Specimen Container	Transport Environment
Preferred	3.6 mL (1.8 mL)	Serum	Clot Activator SST	Room Temperature
Instructions	Recommended for 2 years to 12 years of age. Specimen Type: Red-top tube or gel-barrier tube, Separate Serum from Cells into a screwtop transport container. Stability Requirements: Room temperature 14 days, Refrigerated 14 days, Frozen 3 months. (Freeze/thaw cycles Stable x3)			

GENERAL INFORMATION	
Expected TAT	3-5 Days
CPT Code(s)	86003x8
Service Provided By	 labcorp Oklahoma, Inc.