

**Flow Cytometry on Tissue Biopsy**

Order Name: **FLOW TISS**  
Test Number: 8090065  
Revision Date: 07/06/2022

TEST NAME	METHODOLOGY	LOINC CODE
Flow Cytometry on Tissue Biopsy	Flow cytometry .....	

SPECIMEN REQUIREMENTS				
Specimen	Specimen Volume (min)	Specimen Type	Specimen Container	Transport Environment
Preferred	5x5mm	Tissue	Sterile Container with RPMI Solution	Refrigerated
Alternate 1	5x5mm	Tissue	Sterile Container with Saline	Refrigerated
Instructions	Place a minimum of a 5x5mm tissue sample in a sterile screwtop container with RPMI. Keep specimen Refrigerated. Indicate type of tissue on specimen and requisition. Please deliver to performing laboratory department (flow cytometry) ASAP. <b>DO NOT Place Tissue Into a Fixative Solution.</b> Specimen stability: 48hrs refrigerated			

GENERAL INFORMATION	
Testing Schedule	Tue-Sat
Expected TAT	2 Days
CPT Code(s)	Test Dependent
Service Provided By	 <b>labcorp</b> Oklahoma, Inc.