


**Acid Phosphate Total**

Order Name: **ACID PHS T**  
Test Number: 2000280  
Revision Date: 12/12/2022

TEST NAME	METHODOLOGY	LOINC CODE
Acid Phosphate Total	Enzymatic .....	12173-1

SPECIMEN REQUIREMENTS				
Specimen	Specimen Volume (min)	Specimen Type	Specimen Container	Transport Environment
Preferred	1.5 mL (0.8 mL)	Serum	Clot Activator (Red Top, No-Gel)	Frozen
Instructions	Allow specimen to clot completely. Transfer 1.5mL of serum into an aliquot tube and freeze. <b>CRITICAL FROZEN. Separate specimens must be submitted when multiple tests are ordered.</b> Unacceptable Conditions: Plasma. Non-frozen specimens. Hemolyzed specimens. Stability After separation from cells: Ambient: Unacceptable; Refrigerated: Unacceptable; Frozen: 1 month			

GENERAL INFORMATION	
Expected TAT	5 - 8 days
Clinical Use	This test is used in the evaluation of prostate cancer, Padgett's disease, hyperthyroidism, and Gaucher's disease.
Notes	Labcorp Test Code: 829642
CPT Code(s)	84060
Service Provided By	 <b>labcorp</b> Oklahoma, Inc.