

Beryllium Level

Instructions

Order Name: BERYLLI B
Test Number: 3653400
Revision Date: 12/20/2022

TEST NAME		METHODOLOGY		LOINC CODE
Beryllium Level		Inductively-Coupled Plasma/Mass Spectrometry		metry
SPECIMEN REQ	UIREMENTS			
Specimen	Specimen Volume (min)	Specimen Type	Specimen Container	Transport Environment

Preferred 2.0 mL (0.5 mL) Whole Blood Royal Blue - EDTA (Trace Element Room Temperature Free)

Specimen Type: Royal blue-top (EDTA) metal-free tube

Specimen Storage: SUBMISSION/TRANSPORT (<3 days): Room Temperature. For storage beyond three days, specimen should be refrigerated or

frozen.

Specimen Collection: Use metal-free collection tube.

Special Instructions: Testing referred to MedTox Laboratories Inc MX#60657

GENERAL INFORMATION	
Expected TAT	3-7 days
Notes	Labcorp Test Code: 811752
CPT Code(s)	83018
Service Provided By	labcorp Oklahoma, Inc.