


Methaqualone Quantitative, Urine

Order Name: **Methaqual Ur**
Test Number: 6900779
Revision Date: 04/24/2021

TEST NAME	METHODOLOGY	LOINC CODE
Methaqualone Quantitative, Urine	<u>Gas Chromatography/Mass Spectrometry (GC/MS)</u>	3787-9

SPECIMEN REQUIREMENTS				
Specimen	Specimen Volume (mL)	Specimen Type	Specimen Container	Transport Environment
Preferred	20 mL	Urine, Random	Sterile Screwtop Container	Room Temperature
Instructions	<p>Specimen Type: Use plastic urine drug bottle and tamper-evident seal for forensic specimen. Collection kits are available by request from the laboratory.</p> <p>Specimen Storage: Maintain specimen at room temperature. Stability: If arrival at lab will extend beyond 7 days then refrigerate.</p> <p>Specimen Collection: Urine temperature monitoring is recommended for samples to be tested for medicolegal purposes.</p> <p>Special Instructions: Chain-of-custody documentation is required for samples submitted for pre-employment, random employee testing, and forensic purposes. For other applications, use the standard request form. Please mark chain-of-custody test number on the request form.</p> <p>Specimen Stability: Ambient: Not Available, Refrigerated : Not Available, Frozen: Not Available</p>			

GENERAL INFORMATION	
Notes	Labcorp Test Code: 071332
CPT Code(s)	80375
Service Provided By	 labcorp Oklahoma, Inc.