

**Fanconi Anemia, Group C (FANCC) 2 Mutations**

Order Name: **Fanconi Anem Mut**  
Test Number: 6905009  
Revision Date: 04/15/2017

TEST NAME	METHODOLOGY	LOINC CODE
Fanconi Anemia, Group C (FANCC) 2 Mutations	<u>Polymerase Chain Reaction</u>	

SPECIMEN REQUIREMENTS				
Specimen	Specimen Volume (mL)	Specimen Type	Specimen Container	Transport Environment
Preferred	7 mL (3 mL)	Whole Blood	EDTA (Lavender Top)	Room Temperature
Instructions	<p><b>Notes:</b> 3 mL whole blood, 5 mL amniotic fluid, 10 mg CVS, OR two buccal swabs</p> <p><b>Specimen Type:</b> Lavender-top (EDTA) tube, yellow-top (ACD) tube, sterile plastic conical tube or two confluent T-25 flasks for fetal testing, OR LabCorp buccal swab kit</p> <p><b>Specimen Storage:</b> Maintain specimen at room temperature or refrigerate.</p> <p><b>Specimen Collection:</b> Not Available</p> <p><b>Special Instructions:</b> If cultured cells are needed, an additional 7-12 days may be required. Additional culture fee may be included.</p>			

GENERAL INFORMATION	
Expected TAT	8-15 days
Notes	Labcorp Test Code: 511212
CPT Code(s)	81242
Service Provided By	 <b>labcorp</b> Oklahoma, Inc.