


**Pseudoephedrine Confirmation Urine**

Order Name: **PSEUDO CNF**  
Test Number: 4000138  
Revision Date: 05/14/2020

TEST NAME	METHODOLOGY	LOINC CODE
Pseudoephedrine Confirmation	<a href="#">Liquid Chromatography/Mass Spectroscopy</a>	4004-8
Pseudoephedrine Interpretation	<a href="#">Interpretive information.</a>	

SPECIMEN REQUIREMENTS				
Specimen	Specimen Volume (min)	Specimen Type	Specimen Container	Transport Environment
Preferred	10 mL (1 mL)	Urine, Random	Sterile Container or Monovette	Refrigerated
Instructions	Keep refrigerated. Specimen Stability: Ambient 7 days. Refrigerated 30 days. Frozen 6 months Please provide accurate prescription list to lab.			

GENERAL INFORMATION	
Testing Schedule	Mon - Fri
Expected TAT	3-5 days
Clinical Use	Used to evaluate medication compliance in pain management.
CPT Code(s)	G0480 (alternate:80324)
Service Provided By	 <b>labcorp</b> Oklahoma, Inc.