

Rapid HIV (Exposure)

Order Name: EX RAP HIV

Test Number: 5554025

Revision Date: 05/02/2024

| TEST NAME            | METHODOLOGY       | LOINC CODE |
|----------------------|-------------------|------------|
| Rapid HIV (Exposure) | Rapid Immunoassay |            |

| SPECIMEN REQUIREMENTS |  |               |                     |                       |
|-----------------------|--|---------------|---------------------|-----------------------|
| Specimen              | Specimen Volume (min)  | Specimen Type | Specimen Container  | Transport Environment |
| Preferred             | 3 mL (1 mL)  | Whole Blood   | EDTA (Lavender Top) | Room Temperature      |
| Instructions          | Stability: Whole Blood, Room Temperature 48 hours, Refrigerated 6 days.<br><br>This is part of the Exposure Protocol for non-St. John Health System exposures. Exposure is an unfortunate exchange of blood or body fluid between the "Source" and the "Recipient" by needle stick or other means.<br><br>HIPAA requires the office to give Labcorp Oklahoma, Inc. a phone number and a Name of the person to contact with the stat results. |               |                     |                       |

| GENERAL INFORMATION |   |
|---------------------|---|
| Testing Schedule    | Daily   |
| Expected TAT        | 1 Day   |
| Clinical Use        | HIV Rapid for Exposure  |
| CPT Code(s)         | 86703   |
| Service Provided By |  labcorp<br>Oklahoma, Inc. |