St. John Health System
Lab Catalog

Tapentadol, MS, Ur, MAT

Instructions

Order Name: Tapentadol UR

Test Number: 5195583 Revision Date: 01/09/2024

TEST NAME METHODOLOGY LOINC CODE Tapentadol, MS, Ur, MAT Mass Spectrometry (MS) **SPECIMEN REQUIREMENTS** Specimen Specimen Volume (min) Specimen Type Specimen Container Transport Environment Preferred 30 mL (15 mL) **Sterile Urine Container** Urine, Random **Room Temperature Preservative Free**

Specimen: 30mL(15mL) Random Urine Preservative free Urine Transport tube Room Temperature or Refrigerated.

 $\label{lem:cause} \textbf{Cause for Rejection: Insufficient volume, no ID on container, urine from preservative tube.}$

SUBMISSION/TRANSPORT (<3 days): Room Temperature. For storage beyond three days, specimen should be refrigerated or frozen.

Stability: Room Temperature 5 days, Refrigerated 14 days, Frozen 14 days.

GENERAL INFORMATION	
Expected TAT	5-8 days
CPT Code(s)	80372
Service Provided By	labcorp Oklahoma, Inc.