

Oxycodone Class, MS, Ur, MAT

Order Name: **Oxycodone CI UR**
Test Number: 5195584
Revision Date: 01/09/2024

TEST NAME	METHODOLOGY	LOINC CODE
Oxycodone Class, MS, Ur, MAT	Mass Spectrometry (MS)	

SPECIMEN REQUIREMENTS				
Specimen	Specimen Volume (min)	Specimen Type	Specimen Container	Transport Environment
Preferred	30 mL (15 mL)	Urine, Random	Sterile Urine Container Preservative Free	Room Temperature
Instructions	Specimen: 30mL(15mL) Random Urine Preservative free Urine Transport tube Room Temperature or Refrigerated. Cause for Rejection: Insufficient volume, no ID on container, urine from preservative tube. SUBMISSION/TRANSPORT (<3 days): Room Temperature. For storage beyond three days, specimen should be refrigerated or frozen. Stability: Room Temperature 5 days, Refrigerated 14 days, Frozen 14 days.			

GENERAL INFORMATION	
Expected TAT	5-8 days
CPT Code(s)	80365
Service Provided By	 labcorp Oklahoma, Inc.