

St. John Health System Lab Catalog

Propoxyphene, MS, Ur, MAT

Order Name: **Propoxyphene UR** Test Number: 5195585 Revision Date: 01/09/2024

TEST NAME			METHODOLOGY	LOINC CODE
Propoxyphene, MS, Ur, MAT		Mass Spectrometry (MS)		
SPECIMEN REQU	IREMENTS			
Specimen	Specimen Volume (min)	Specimen Type	Specimen Container	Transport Environment
Preferred	30 mL (15 mL)	Urine, Random	Sterile Urine Container Preservative Free	Room Temperature
Instructions	Specimen: 30mL(15mL) Random Urine Preservative free Urine Transport tube Room Temperature or Refrigerated. Cause for Rejection: Insufficient volume, no ID on container, urine from preservative tube. SUBMISSION/TRANSPORT (<3 days): Room Temperature. For storage beyond three days, specimen should be refrigerated or frozen. Stability: Room Temperature 5 days, Refrigerated 14 days, Frozen 14 days.			

GENERAL INFORMATION	
Expected TAT	5-8 days
CPT Code(s)	80367
Service Provided By	Oklahoma, Inc.