

Specimen Volume (min)

Specimen

Instructions

Order Name: Buprenor Confirm UR

Transport Environment

Test Number: 5195594
Revision Date: 01/09/2024

Buprenorphine Confirmation, Ur

Liquid Chromatography/Tandem Mass Spectrometry

SPECIMEN REQUIREMENTS

METHODOLOGY

Liquid Chromatography/Tandem Mass Spectrometry

Preferred 30 mL (15 mL) Urine, Random Sterile Urine Container Room Temperature
Preservative Free

Specimen: 30mL(15mL) Random Urine Preservative free Urine Transport tube Room Temperature or Refrigerated.

Cause for Rejection: Insufficient volume, no ID on container, urine from preservative tube.

Specimen Type

SUBMISSION/TRANSPORT (<3 days): Room Temperature. For storage beyond three days, specimen should be refrigerated or frozen.

Specimen Container

Stability: Room Temperature 14 days, Refrigerated 14 days, Frozen 6 weeks.

GENERAL INFORMATION	
Expected TAT	5-8 days
CPT Code(s)	80348
Service Provided By	labcorp Oklahoma, Inc.