

Order Name: Antidepre Confirm UR Test Number: 5195617 Revision Date: 01/11/2024

		METHODOLOGY	LOINC CODE
Antidepressants Confirm, UR		Liquid Chromatography/Tandem Mass Spectrometry	
REMENTS			
Specimen Volume (min)	Specimen Type	Specimen Container	Transport Environment
30 mL (15 mL)	Urine, Random	Sterile Urine Container Preservative Free	Room Temperature
Cause for Rejection: Insufficien SUBMISSION/TRANSPORT (ent volume, no ID on container, ur (<3 days): Room Temperature. Fo	rine from preservative tube. or storage beyond three days, specimen sho	•
	REMENTS Specimen Volume (min) 30 mL (15 mL) Specimen: 30mL(15mL) Rande Cause for Rejection: Insufficien SUBMISSION/TRANSPORT (nfirm, UR REMENTS Specimen Volume (min) Specimen Type 30 mL (15 mL) Urine, Random Specimen: 30mL(15mL) Random Urine Preservative free Urine Cause for Rejection: Insufficient volume, no ID on container, ur SUBMISSION/TRANSPORT (<3 days): Room Temperature. For	REMENTS Specimen Volume (min) Specimen Type 30 mL (15 mL) Urine, Random Sterile Urine Container

GENERAL INFORMATION	
Expected TAT	5-8 days
CPT Code(s)	80338, 80334, 80337
Service Provided By	Oklahoma, Inc.