


Vancomycin Trough

Order Name: **VANCO TROU**
Test Number: 4005950
Revision Date: 06/04/2021

TEST NAME	METHODOLOGY	LOINC CODE
Vancomycin Trough	Enzyme Multiplied Immunoassay Technique	4092-3

SPECIMEN REQUIREMENTS				
Specimen	Specimen Volume (min)	Specimen Type	Specimen Container	Transport Environment
Preferred	1 mL (0.5)	Plasma	Lithium Heparin PST (Light Green Top)	Refrigerated or Frozen
Alternate 1	1 mL (0.5)	Serum	Clot Activator SST	Refrigerated or Frozen
Instructions	Trough: Draw specimen immediately preceding next dose. Separate plasma/serum from cells/gel immediately. Stability: Refrigerated 7 days, Frozen 14 days.			

GENERAL INFORMATION	
Testing Schedule	Daily
Expected TAT	1-2 days
Clinical Use	Useful for optimizing drug dosage and assessing toxicity.
CPT Code(s)	80202
Service Provided By	 labcorp Oklahoma, Inc.