

Rotavirus Screen

Order Name: C ROTA SC

Test Number: 6100250

Revision Date: 09/06/2024

TEST NAME	METHODOLOGY	LOINC CODE
Rotavirus Screen	Enzyme Immunoassay	

SPECIMEN REQUIREMENTS				
Specimen	Specimen Volume (mL)	Specimen Type	Specimen Container	Transport Environment
Preferred	3mL	Stool, Raw	Sterile Screwtop Container	Refrigerated
Instructions	Specimen is good for 72 hours refrigerated.			

GENERAL INFORMATION	
Testing Schedule	Daily
Expected TAT	1 Day
Clinical Use	Detects rotavirus in stool
CPT Code(s)	87425
Service Provided By	<div>labcorp Oklahoma, Inc.</div>