

Nut Allergen Panel

Order Name: **A NUTS PNL**
Test Number: **5616500**
Revision Date: **08/28/2023**

TEST NAME	METHODOLOGY	LOINC CODE
Peanut Allergen	ImmunoCAP	15917-8
Walnut Allergen	ImmunoCAP	16074-7
Pecan Allergen (Food)	ImmunoCAP	15773-5
Almond Allergen	ImmunoCAP	15527-5
Hazelnut Allergen	ImmunoCAP	15766-9
Brazil Nut Allergen	ImmunoCAP	15586-1
Cashew Allergen	ImmunoCAP	15607-5
Pistachio Allergen	ImmunoCAP	15949-1

SPECIMEN REQUIREMENTS				
Specimen	Specimen Volume (min)	Specimen Type	Specimen Container	Transport Environment
Preferred	3.2 mL (1.6 mL)	Serum	Clot Activator SST	Room Temperature
Instructions	Specimen Type: Red-top tube or gel-barrier tube, Separate Serum from Cells into a screwtop transport container. Stability Requirements: Room temperature 14 days, Refrigerated 14 days, Frozen 3 months. (Freeze/thaw cycles Stable x3)			

GENERAL INFORMATION	
Expected TAT	3-5 Days
CPT Code(s)	86003x8
Service Provided By	 labcorp Oklahoma, Inc.