

**Meconium Drug Screen (10 Panel)**

Order Name: **DRUG F S10**  
Test Number: 4066450  
Revision Date: 10/01/2022

TEST NAME	METHODOLOGY	LOINC CODE
Meconium Drug Screen (10 Panel)		

SPECIMEN REQUIREMENTS				
Specimen	Specimen Volume (min)	Specimen Type	Specimen Container	Transport Environment
Preferred	5g (1g)	See Instructions	See Special Instructions	Refrigerated
Instructions	<p><b>CHAIN OF CUSTODY TEST !</b> - Please contact RML Client Services to arrange for collection supplies and instruction before submission of specimen for testing.</p> <p>5g(1g) Meconium in an appropriate Chain of Custody packet, Refrigerated.</p> <p><b>Specimen Type:</b> Meconium Collection Kit (PeopleSoft#96663) OR clean, leakproof polypropylene container.</p> <p><b>Specimen Storage:</b> SUBMISSION/TRANSPORT (&lt;3 days): Room Temperature. For storage beyond three days, specimen should be refrigerated or frozen.</p> <p><b>Specimen Collection:</b> Specimens from different voidings may be pooled if necessary</p> <p><b>Special Instructions:</b> Testing referred to MedTox Laboratories Inc MX#4051</p>			

GENERAL INFORMATION	
Expected TAT	3-6 Days
Notes	Labcorp Test Code: 809313
CPT Code(s)	80307
Service Provided By	 <b>labcorp</b> Oklahoma, Inc.