

Mold Allergy Panel

Order Name: **AO MOLD PN**
Test Number: 2929845
Revision Date: 05/06/2024

TEST NAME	METHODOLOGY	LOINC CODE
<u>Alternaria alternata Allergen</u>	<u>ImmunoCAP</u>	15530-9
<u>Epicoccum purpurascens Allergen</u>	<u>ImmunoCAP</u>	15700-8
<u>Aspergillus fumigatus Allergen</u>	<u>ImmunoCAP</u>	15549-9
<u>Cladosporium herbarum Allergen</u>	<u>ImmunoCAP</u>	6075-6
<u>Fusarium moniliforme Allergen</u>	<u>ImmunoCAP</u>	15722-2
<u>Aureobasidium pullulans Allergen</u>	<u>ImmunoCAP</u>	15964-0
<u>Penicillin Chrysogonium Allergen</u>	<u>ImmunoCAP</u>	15924-4

SPECIMEN REQUIREMENTS				
Specimen	Specimen Volume (min)	Specimen Type	Specimen Container	Transport Environment
Preferred	1mL (0.5mL)	Serum	Clot Activator SST	Refrigerated
Instructions	Specimen Stability: Room Temperature: Unacceptable, Refrigerated: 7 days, Frozen 30 Days			

GENERAL INFORMATION	
Testing Schedule	Tue-Sat
Expected TAT	2-4 Days
CPT Code(s)	86003x7
Service Provided By	 labcorp Oklahoma, Inc.