

Testosterone Total, LCMS

Order Name: **TEST T LCMS**
Test Number: 3602651
Revision Date: 06/28/2023

TEST NAME	METHODOLOGY	LOINC CODE
Testosterone, Total	<u>Liquid Chromatography/Tandem Mass Spectrometry</u>	2986-8

SPECIMEN REQUIREMENTS				
Specimen	Specimen Volume (mL)	Specimen Type	Specimen Container	Transport Environment
Preferred	1 mL (0.5 mL)	Serum	Clot Activator SST	Refrigerated
Instructions	<p>Recommended for Women, children under 13 years of age, hypogonadal males, or patients undergoing anti-androgenic therapies.</p> <p>Specimen: Serum - Red-top tube or serum separator (SST)</p> <p>Collection: Serum should be transferred from cells within one hour of collection and transferred to a plastic transport tube.</p> <p>Preferred Specimen Volume: 1.0 mL (0.2 mL min - does not allow for repeat testing)</p> <p>Storage/ Transport Temperature: Room temperature, refrigerated or frozen</p> <p>Stability: Room Temp: 8 days; Refrigerated: 7 days; Frozen: 5.58 years</p> <p>Cause for Rejection: Improperly labeled specimens, specimens collected tube other than red-top or SST.</p>			

GENERAL INFORMATION	
Testing Schedule	Mon-Fri
Expected TAT	2-5 Days
CPT Code(s)	84403
Service Provided By	 labcorp Oklahoma, Inc.